AUTOCLAVE A VAPORE
STEAM AUTOCLAVE
AUTOCLAVE A' VAPEUR

STERIL MAXIMUM

CE 0068

RUA DA LAMEIRA, 115
4445-490 ERMSINDE
912313018
geral@portodonto.com
www.portodonto.com

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
1.3 INFORMATION TABLE (Annex D - EN13060)


REQUIREMENTS	B	S
Dynamic sterilizer chamber pressure	X	X
Air outlet	X	X
Empty chamber	X	X
Solid load	X	X
Small porous objects	X	
Small porous loads	X	
Full porous load	X	
B type hollow load	X	
A type hollow load	X	
Multiple wrap	X	
Drying, solid load	X	X
Drying, porous load	X	X
Residual air		

1.4 GRAPHIC EXAMPLES OF VARIOUS TYPES OF CYCLE

type B cycle

type S cycle

1.1 INTENDED USE

Dear Customer,

The sterilizer is a device designed for water steam sterilization of small sized tools and professional equipment, and is widely used for medical purposes by general practitioners and dentists, in salons that deal with personal hygiene and body care and also in veterinary surgeries. It is also used for sterilizing materials and equipment that come into contact with blood or physiological liquids, such as instruments used by beauticians, tattoo artists, piercers and hairdressers. The very specific sterilization loads used in these sectors of applications require different performance features for the sterilization cycles. It is, therefore, advisable that this device is used by professional personnel duly trained.

It is of fundamental importance for the sterilizer and the respective equipment to be used solely for the type of product for which they were designed. We therefore request you to consult the Declaration of Conformity of this appliance: the Class to which the appliance belongs is indicated in the "Category" box. The "Sterilization Table" (Chap. 11) gives all the information necessary for establishing the type of cycle to be used for sterilizing the various instruments.

Misuse of or failure to follow the instructions in this manual can cause danger to the user.

1.2 GENERAL WARNINGS

- We recommend that you read the Instruction Manual carefully before using the device, to ensure that the operations required are carried out correctly: **DO NOT** carry out operations other than those described in this booklet. The Manufacturer declines all responsibility for direct or indirect damage to objects, persons or animals deriving from improper use of the appliance.
- The machine must be used by responsible adults only.
- Position the machine in a place where it is not accessible to children.
- Install the machine in a position where it is easy to reach the plug.
- Do not use the machine near flammable or explosive sources.
- Use the machine in dry protected areas.
- Check the condition of the power cable periodically: do not use the appliance if the power cable is not perfectly intact.
- Do not carry out maintenance with the machine running or if it is plugged into the power socket.
- Do not approach the machine with flammable material.
- Always use personal protection devices, in compliance with the applicable directives.
- Do not use the appliance for purposes other than those mentioned in this Instruction Manual.
- Read the paragraph concerning the technical features carefully before starting up the appliance.
- For your safety, make sure you pay great attention to the instructions given below.
- The sterilizer must be used in a controlled environment and by trained personnel.
- Do not perform maintenance or cleaning operations while the machine is operating or plugged in.
- Do not stand in front of the sterilizer when in use.
- The sterilizer has not been designed to sterilize food products. Do not insert food.

IT IS STRICTLY FORBIDDEN TO USE THIS STERILIZER UNDER THE INFLUENCE OF MEDICINES, DRUGS OR ALCOHOL

1. TABLE OF CONTENTS

1.	TABLE OF CONTENTS
1.1	INTENDED USE
1.2	GENERAL WARNINGS
1.3	INFORMATION TABLE (Annex D - EN13060)
1.4	GRAPHIC EXAMPLE OF VARIOUS TYPES OF CYCLE
2.	SAFETY
2.1	SAFETY SYMBOLS
2.2	SAFETY DEVICES
3.	PACKAGING, STORAGE AND DISPOSAL
3.1	DISPOSAL AND/OR SCRAPPING
4.	FIRST INSTALLATION
4.1	WATER CONNECTIONS
5.	ACCESSORIES
6.	CONTROL PANEL
6.1	DISPLAY
6.2	KEYBOARD
6.3	INDICATOR ICONS
7.	PRINTER
7.1	CHANGING THE PRINTER PAPER ROLL
7.2	PRINT HEAD CLEANING
7.3	PRINTOUTS
8.	USER MENU
8.1	ENGLISH
8.2	TIME (hh:mm:ss)
8.3	DATE (dd:mm:yy)
8.4	ADD WATER...
8.5	CYCLES RE-PRINT
8.6	CYCLES RE-PRINT AL0000
8.7	EMERGENCY ON/OFF
9.	FIRST TIME START-UP
9.1	DOOR CLOSURE
9.2	TEST CYCLE
10.	WATER LOADING AND DISCHARGE
10.1	CLEAN WATER: TANK FILLING VIA PUMP
10.2	CLEAN WATER: AUTOMATIC TANK FILLING
10.3	DIRTY WATER: MANUAL TANK DRAINING
10.4	DIRTY WATER: AUTOMATIC TANK DRAINING
10.5	CLEAN WATER: TANK DRAINING AND CLEANING
10.6	WATER QUALITY TABLE
11.	STERILIZATION TABLE
11.1	AUTOCLAVE B6 TYPE
11.2	AUTOCLAVE B18 TYPE
11.3	AUTOCLAVE B23 TYPE
11.4	NIGHT CYCLE
12.	TEST CYCLES
12.1	"BOWIE & DICK" CYCLE
12.2	"HELIX TEST" CYCLE
12.3	VACUUM TEST CYCLE
12.4	BIOLOGICAL TEST
13.	RECOMMENDATIONS FOR STERILIZATION
14.	ALARMS AND ERRORS
14.1	ALARMS
14.2	ERRORS
15.	MAINTENANCE
15.1	DAILY MAINTENANCE
15.2	WEEKLY MAINTENANCE
15.3	QUARTERLY MAINTENANCE
15.4	ANNUAL MAINTENANCE/VALIDATION
16.	TECHNICAL FEATURES
17.	WARRANTY

2. SAFETY

2.1 SAFETY SYMBOLS

	ATTENTION: READ THE USER MANUAL INSTRUCTIONS CAREFULLY
	VOLTAGE
	ATTENTION HIGH TEMPERATURE
	EARTH CONNECTION

2.2 SAFETY DEVICES

This equipment includes the following safety devices:

1. Three micro-switches controlling the door and the automatic locking system: they are independent of each other and ensure that the door is closed and properly locked. If problems arise, the user is warned by an alarm and the next cycle cannot start. If any problem arises while the cycle is running, the microprocessor stops the process and immediately reduces the machine pressure level.
2. Two different mechanical thermostats ensure that the temperature of the various components does not exceed the pre-set value. Thermostats are manually reset.
3. Four electronic temperature sensors continually monitor all critical points on the machine, to prevent over temperature errors during the operating process.
4. A pressure relief safety valve protects against the risk of explosion. Safety valve inspection is reserved for Regulatory authorities and/or authorised technicians and is governed by specific law, in force in the country of installation. Please refer to the user manual of the valve provided.
5. An electronic pressure transducer checks all the solenoid valves and opens them if an overpressure is detected.

3. PACKAGING, STORAGE AND DISPOSAL

The cardboard packaging used for transporting the sterilizer IS NOT STERILE.

The sterilizer is fragile and should be transported with extreme care.
DO NOT TURN UPSIDE DOWN.

The sterilizer is packed with all accessories placed inside the chamber. It is wrapped in a protective polyethylene bag and then placed inside a cardboard box. To protect it against accidental impact, it is also padded with polystyrene, cardboard, or polyurethane.

Upon delivery, check the conditions of the packaging box and the sterilizer. In case of damage, contact your dealer, the shipper and the carrier who performed transport.

The user is recommended to retain the packaging for the period of the warranty: if the equipment is returned for repair without the original packaging, the user will be charged for new packaging at time of reshipment.

Unpacking and locating the sterilizer:

At least two people are required to remove the equipment from the box. Follow the instructions given below:

- Open the box and remove the staples to avoid getting scratched or cut while removing the equipment.
- Carefully read the instructions for use.
- Remove the machine from its box by taking hold of the straps on the sides. This will prevent putting strain on plastic parts.
- Place the machine on a perfectly horizontal level surface with a load capacity of at least 70 kg.
- Connect the plug to a Schuko plug with safety ground.
Do not replace the original plug with other types.
Do not make additional connections.
Do not connect to multiple plugs. In addition, ensure that the installation to which the sterilizer is connected complies with all relevant legislation and can sustain the specified load (paragraph 16).
- Open the door by pushing the DOOR button.
- Take the accessory kit out of the box and switch the equipment off.
- CAREFULLY READ PARAGRAPHS 4 AND 4.1 BEFORE STARTING NORMAL OPERATION.
- Keep in a dry and protected area at a temperature of 5-30°C.

3.1 DISPOSAL AND/OR SCRAPPING

Consult local codes and ordinances for proper disposal and/or scrapping of any component and other consumable goods (packaging, water, etc.).

Directive 2002/96/EC (Waste Electrical and Electronic Equipment - WEEE): information for users. This product complies with EU Directive 2002/96/ EC. The crossed-out wastebasket symbol on the appliance means that at the end of its useful lifespan, the product must be disposed of separately from ordinary household wastes. The user is responsible for delivering the appliance to an appropriate collection facility at the end of its useful lifespan. Appropriate separate collection to permit recycling, treatment and environmentally compatible disposal helps prevent negative impact on the environment and human health and promotes recycling of the materials making up the product. For more information on available collection facilities, contact your local waste collection service or the distributor where you bought this appliance.

4. FIRST INSTALLATION



Correct installation of the sterilizer is key to ensuring its proper operation. Below is an installation checklist:

1. The device must be installed in an area accessible only to authorised personnel.
2. The working environment must be properly lit and adequately ventilated.
3. Position the equipment on a perfectly flat and horizontal surface with a minimum load capacity of 70 kg. The sterilizer is pre-levelled when shipped. The sterilization chamber is slightly inclined to the rear of the equipment. Leave at least 5 cm between the wall and the back of the sterilizer.
4. Position the sterilizer so as to leave sufficient room for inspecting and cleaning the sterilization chamber.
5. Do not install the sterilizer close to washbasins or taps: the cover of this equipment is not watertight.
6. Do not install the equipment close to heat sources (i.e. other sterilizers, ovens or similar).
7. In order to avoid any damage to people, objects or animals, the machine must be put in a position allowing the flow coming from the safety valve to flow in a safe place.

4.1 WATER CONNECTIONS

The equipment water connections are of primary importance for its correct operation. Always leave some space free (at least 5 cm) all around the sterilizer, especially from the rear and the wall.

Follow instructions as detailed in the photo.

ATTENTION:

If the automatic water filling system "D" is used, you MUST connect the condensate drainage pipe "B" to a drain line and not to the condensate container.

5. ACCESSORIES


The tray holder is supplied together with four trays and the following accessories.

	CLEANING SPONGE The sponge should be used as indicated in paragraph 15.
	SILICON-INSULATED GLOVE Use the glove to extract the hot trays. or TOOL FOR TRAY EXTRACTION Use the right side to extract the hot trays.
	1 FILL HOSE Install the filter in the water container and secure the other end to the hose connection at the top, left-hand side of the chamber rim: press the PUMP/ WATER button to start filling.
	2 DRAIN TUBES: Use one tube as described in paragraph 10.3, to drain the used water. The second tube should be connected to the hose connection "B", as described in paragraph 4.1.

6. CONTROL PANEL

In order to understand and use the controls correctly, please refer to paragraphs 6.1, 6.2 and 6.3.

6.1 DISPLAY

6.2 KEYBOARD

DOOR	This button opens the door at the end of the cycle. By keeping this pressed you can see the total STR+DRY time and the cycle count.
SELECT	This button selects the desired cycle. Also used to enter the USER MENU (paragraph 8).
PUMP / WATER	Used to fill water in the tank. Standard autoclaves: hold it down for at least 5 sec. for the resin reset procedure (Chapter 14.2 Reset Resins) Autoclaves with water quality detection system: if pressed after a "BAD WATER" message, reactivates automatic tank filling (point 10.2).
START/STOP	Starts the cycle. While the cycle is running, holding the button down for longer than 1 sec. will stop the cycle.

SELECT + START	Hold it down, for at least 10 sec., for the alarm reset (Paragraph 14)
START 10 sec.	Hold it down, with the door open and for at least 10 sec., for the BACT. FILTER reset (Paragraph 14.2)

6.3 INDICATOR ICONS

	This icon is displayed shortly after a cycle has started and indicates that the machine door is locked.
	This icon indicates that a cycle is running. The black triangle rotates in the circle.
	This icon warns that printer is out of paper, or the printer cover is not properly closed. The machine can execute its cycle without any risk.
	Type of cycle selected: B, S or N type (see Table on paragraph 1.3).
	The sterilizers are equipped with software which is able to communicate and to write data onto a USB key. The sterilizers have a USB connector already installed on the machine. When the USB key is inserted, the USB icon appears on the display. It is always necessary to read the specific enclosed manual carefully to prevent data saving errors. <u>DATA SAVING BEGINS AND ENDS ONLY AND EXCLUSIVELY IF THE USB KEY IS ALREADY INSERTED IN THE CONNECTOR BEFORE THE CYCLE IS STARTED.</u>
	This icon warns the user that clean water is at its minimum level: it is not possible to start a new cycle. Fill the machine with water before starting the next cycle (see paragraphs 10.1 and 10.2).
	The clean water tank is full. If the PUMP button is pressed again, the display will show that it is not possible to continue filling the tank with water.
	While this icon is displayed, no cycle can start: the dirty water tank should be completely drained. <u>Allow all water to drain out before closing the draining tap.</u> (see paragraphs 10.3 and 10.4).

7. PRINTER





On starting any cycle, the sterilizer starts printing all values for the selected cycle, plus the machine model and serial number (paragraph 7.3). When the sterilizer cycle is completed, the printer stops printing: cut off the printout using the built-in cutter (simply pull the paper upward). If the cover is not properly closed, or if the printer is out of paper, the relevant indicator icon will be displayed to warn the user (paragraph 6.3).

Refer to the instructions provided by the local health authority for printout-filing.

For a correct and long-lasting storage of the paper printouts, it is necessary to keep them away from light and heat sources.


7.1 CHANGING THE PRINTER PAPER ROLL

To change the printer paper roll: open the cover, insert a new thermal paper roll (the maximum roll diameter accepted is 57mm) in the space provided and feed the paper over the cover roll while closing the cover. Use only thermal paper. Make sure the thermal paper is correctly oriented: if the paper is inserted the wrong way up, the printouts will be blank.

7.2 PRINT HEAD CLEANING

If the printouts are difficult to read, the print head should be cleaned using a lint free cloth moistened with alcohol: open the printer cover, remove the paper roll and clean the head (see the picture). Use compressed air to clear dust from the interior of the printer.

7.3 PRINTOUTS

How to read and interpret the information provided in the printouts.

Model.	STAR	Software release.
Serial number.	RELEASE: 0.02 S/N: STAR1234	Cycle start date.
Selected cycle.	DATA: 28:12:04 C/N: 0312	Incremental number of cycles started.
Column showing time: hour, minutes, seconds.	CICLO: BOWIE & DICK 134° 2.06bar 3.30min	Info on the cycle underway: minimum cycle temperature and pressure values.
Column showing temperature values in degrees centigrade (°C).	START hh:mm:ss °C bar	Pressure column.
Current cycle phase.	12:42:25 015.7 0.019 12:47:50 025.3 -0.801 12:55:41 099.7 0.202 12:58:20 068.4 -0.701 13:01:45 103.5 0.200 13:06:32 068.7 -0.701	Cycle result.
Operator's signature.	STERILIZZAZIONE 13:14:54 135.0 2.207 13:15:53 134.9 2.200 13:16:52 134.6 2.181 13:17:52 134.6 2.179 ASCIUGATURA 13:18:53 134.2 2.171 13:19:51 108.8 0.136 13:20:51 095.9 -0.383 13:21:50 090.4 -0.126 FINE CICLO: OK OPERATORE:	

8. USER MENU

To access the user menu, proceed as described below:

- Switch OFF the sterilizer.
- Press and hold SELECT then switch ON the sterilizer: release the SELECT button only when the preset language (e.g.: ENGLISH) appears on the display.

Within the user menu, the pushbuttons are used in the following manner:

- START = proceed to the next page.
- DOOR and PUMP = used to modify the settings in the pages.
- SELECT = when pressed until the confirmation sound is given, it saves the pre-selected value (not all the pages need to save the pre-selected value).

8.1 ENGLISH

Use DOOR and PUMP buttons to select the required language. Press the START/STOP button to continue...

8.2 TIME (hh:mm:ss)

This page allows you to set the time. Press the DOOR button to increase the value and the PUMP button to decrease it; use the SELECT button to move the cursor to the value to be changed. Press the START/STOP button to continue...

8.3 DATE (dd:mm:yy)

This page allows you to set the date. Press the DOOR button to increase the value and the PUMP button to decrease it; use the SELECT button to move the cursor to the value to be changed. Press the START/STOP button to continue...

8.4 CHARGE WATER BY...

This page is used to select the water filling mode (...USING PUMP or ...MANUALLY or ...FROM OSMOSIS or ...FROM DEMINERALIZER). Read carefully the paragraph 10.2.
Press the START/STOP button to continue...

8.5 CYCLES RE-PRINT

This function allows you to save the last cycles executed by the machine with a single action. These cycles can be either printed from the internal or external printer by holding down the SELECT button for 5 sec. Press the START/STOP button to continue...

8.6 CYCLES RE-PRINT AL0000

This function allows you to save the last cycles not completed. These cycles can be either printed from the internal or external printer by holding down the SELECT button for 5 sec. Press the START/STOP button to continue...

8.7 EMERGENCY ON/OFF

The "Emergency ON" setting is reserved for specialised technical personnel.

Press the START/STOP button to return to the language setting page; hold it down to quit the user menu.
You can exit the USER MENU at any time by holding down the START/STOP button.

ENGLISH

9. FIRST TIME START-UP

USING THE SYSTEM FOR THE FIRST TIME

After installing the sterilizer (paragraphs 4 and 4.1), make sure that all water connections are watertight. Proceed as detailed below:

1. Start the equipment using the main switch.
2. Enter the USER MENU (paragraph 8) and configure the system in accordance with the selected filling mode (Refer to paragraph 4.1).
3. Fill the tanks with water: once the clean water MIN. LEVEL icon disappear, the test cycle can start.
4. Insert at least the tray holder and start the cycle by using the START button: the sterilizer will run automatically until the end of the cycle, a warning sound will be heard and the display will show "END CYCLE" and "OPEN THE DOOR".

9.1 DOOR CLOSURE

The procedure to close/open the door is described below:

Procedure to close: Keep the door pushed in until the display shows the door locked icon (paragraph 6.3). In case of closure with handle, lift the handle, close the door and release the handle when the door is closed. To close the door, it is NOT necessary to press the START button as the motor starts automatically.

To start the cycle, press the START button once the door is closed.

Procedure to open: To open the door, it is necessary to press the DOOR button. To guarantee maximum operational safety, the sterilizer will prevent the door from opening until all residual pressure inside the sterilizer is discharged. In this case, simply press the DOOR button again.

9.2 TEST CYCLE


The purpose of the test cycle is to ensure that the sterilizer is in working order, has not suffered any damage during transport and is not subject to any operational faults due to technical causes. It is recommended that this test is based on 134°C cycle.

1. Close the door as described in paragraph 9.1 then, with the door closed, press the START button.
2. Cycle 134°C:
 - While the cycle is running, the display will show the temperature – pressure– time remaining– selected cycle temperature – type of cycle - current cycle phase - any warning icons.
 - Insert the tray holder and the trays inside the machine as indicated at paragraph 12.1. The vacuum pump starts removing air from the chamber until the pre-set value is reached (VACUUM 1), after which the rise phase will start (PREHEATING). During this phase a slight buzz coming from the water injection pump (used to inject steam into the chamber) can be heard.
 - Next the pre-set vacuum degree is reached the EXHAUST phase will start; another two vacuum phases and two preheating phases will be carried out until the cycle pre-set values are reached (paragraph 11).
 - The STERILIZATION phase will now start. During the exposure period, pressure and temperature will be continually software-controlled to ensure that the sterilisation is successful. Alarms will trip in the event of problems (see paragraphs 14 and 14.1).
 - On completion of the sterilization phase, the DRYING phase will begin: pressure inside the chamber will be released and the final vacuum phase will start (the vacuum pump will evacuate any steam to considerably improve the final drying of the sterilized instruments).
 - The door can only be opened when the display shows "END CYCLE" and "OPEN THE DOOR".
 - In the event of residual positive or negative pressure, the door lock system will not open (proceed as indicated in paragraph 6.2).
 - Remove the instruments by using the special handle supplied with the equipment: **use protective gloves to avoid getting burnt.**

10. WATER TANKS FILLING AND DRAINING

10.1 CLEAN WATER: TANK FILLING VIA PUMP

Total clean water tank capacity is 4.5 litres. The autoclave is supplied pre-programmed to perform tank filling via internal pump. To fill the tank, connect the supplied hose (point 5) to the hose connection at the top left of the boiler and press the PUMP button: the charging pump will activate for a max time of 380".

After 380" or when the corresponding icon reaches maximum level (point 6.3), the display will light up and the pump will stop automatically. It is possible to stop the filling operation at any time by simply pressing the PUMP button. Cycle start-up will be inhibited in the event of minimum level as signaled on the display by the corresponding icon.

For autoclaves with water quality detection system (OPTIONAL): during filling, water quality is monitored and various messages can appear on the display.

Autoclaves with water quality detection system (OPTIONAL) TANK FILLING VIA PUMP		
MESSAGE	MEANING	SOLUTION
"POOR"	Means that the quality of the water being filled in the autoclave is at the limits of acceptability. See chapt. 10.6.	Acquire water with better quantity than the currently used water <u>as soon as possible</u> to avoid damage to the autoclave.
"BAD"	The pump stops and tank filling is interrupted, even if the MIN. LEV. icon is still present. Means that the quality of the water being filled in the autoclave is outside of the limits of acceptability. See chapt. 10.6.	Using water different from that previously used, press the PUMP button again: After 15 seconds, if water is still not good quality, the pump will stop once again. After a couple attempts, you <u>must</u> replace water being used with better quality water: continued use of unsuitable water will void the warranty on the autoclave.
NOTE:	The autoclave preventively records water quality. The use of unsuitable quality water shall void the product warranty.	

10.2 CLEAN WATER: AUTOMATIC TANK FILLING

If water loading is automatic and occurs by means of the mains water supply and osmosis system, the water supply system must be equipped with a reverse flow prevention device in compliance with IEC 61770.

If a water purification system (Demineralizer or OSMOSIS system) is present, the water connections as described in paragraph 4.1 must be used and the software must be programmed with the pre-selected setting. Proceed as described below:

1. Switch OFF the sterilizer.
2. Keep SELECT pressed and switch ON the sterilizer: release the SELECT button only when the pre-set language (e.g. ENGLISH) appears on the display.
3. Press START to reach the page in which it is possible to select the clean water tank filling mode. The various possibilities are:
 - **ADD WATER...USING PUMP** (the sterilizer takes in the water from the front: paragraph 10.1).
 - **ADD WATER...MANUALLY** (ATTENTION: for technical use only!)
 - **ADD WATER...FROM OSMOSIS** (the machine needs an osmotic purification system).
 - **ADD WATER...FROM DEMINERALIZER** (the sterilizer requires resin-based purification system).
4. Set the mode by pressing PUMP, then keep SELECT pressed for 5 seconds: the long sound confirms the value is saved in the software. Switch OFF the sterilizer.

5. Read and comply with the instructions of the purification system you are installing.

Once the installation of the purification system is completed, turn the machine on and verify all the water connections while the tank is being filled.


THE MANUFACTURER DECLINES ALL RESPONSIBILITY FOR PROBLEMS DERIVING FROM INSTALLATION AND CONNECTIONS CARRIED OUT BY UNAUTHORIZED PERSONNEL.

Autoclaves with water quality detection system (OPTIONAL) AUTOMATIC TANK FILLING		
MESSAGE	MEANING	SOLUTION
"POOR"	Means that the quality of the water entering into the autoclave tank is at the limits of acceptability. See chapt. 10.6.	Acquire water purification system spare parts and replace them <u>as soon as possible</u> to prevent damage to the autoclave (operation intended for a technician).
"BAD"	The autoclave has stopped tank filling (even if the MIN. LEV. icon is still present). Means that the quality of the water entering into the autoclave tank is outside the limits of acceptability. See chapt. 10.6.	Press the PUMP button to reactivate automatic filling: After 15 seconds, if water is still not good quality, the autoclave will stop filling once again. After a couple attempts, you <u>must</u> intervene on the purification system, replacing the necessary filters (operation intended for a technician): continued use of unsuitable water will void the warranty on the autoclave.
NOTE:	The autoclave preventively records water quality. The use of unsuitable quality water shall void the product warranty.	

10.3 DIRTY WATER: MANUAL TANK DRAINING

The capacity of the dirty water tank is 4.5 litres. To empty the tank, use the silicon tube supplied with the machine, by connecting it to the discharge joint (see photo). Push the joint towards the sterilizer and wait for all the water to drain out before disconnecting the tube. Pull the joint and disconnect the tube. Dirty water may contain contaminated residues and it is therefore recommended that protective gloves are used while draining water from the equipment.

Water discharge joint: use the silicon tube supplied with the machine by connecting it into the discharge joint. Push the joint towards the autoclave and wait for all the water to drain before disconnecting the tube.

REMOVE TUBE AFTER DRAINING.

10.4 DIRTY WATER: AUTOMATIC TANK DRAINING

To eliminate the need to empty the dirty water tank manually, use the tap at the rear of the sterilizer. Connect a hose to the tap, secure it with a metal hose clamp, and then insert the other end of the hose into a container or connect it directly to a water drain (paragraph 4.1). Leave the tap open: the dirty water will be drained automatically into the collection container or into the water drain.

10.5 CLEAN WATER: TANK DRAINING AND CLEANING

In order to empty the tank and gain access for cleaning purposes, proceed as indicated in paragraph 4.1. This operation should be performed only after disconnecting the water supply line from the mains (if the sterilizer is connected to an external water purifier).

10.6 WATER QUALITY TABLE (UNI EN 13060)

The water to be used must be purified or distilled and must meet the requirements specified in the table. Do not use isopropyl alcohol or unsuitable liquids. Do not add any additives or chemicals to water. Do not reuse the drain water from the tank.

UNI EN 13060 – annex C	Maximum value
Evaporation residue	10 mg/l
Silicon oxide (SiO ₂)	1 mg/l
Iron	0.2 mg/l
Cadmium	0.005 mg/l
Lead	0.05 mg/l
Heavy metal residues (except iron, cadmium and lead)	0.1 mg/l
Chloride	2 mg/l
Phosphate	0.5 mg/l
Conductivity (at 20°C)	15 µs/cm
pH value	5 to 7
Look	colourless, clean, without sediment
Hardness	0.02 mmol/l

NOTE: The use of water containing contaminants at concentrations exceeding those given in the table can considerably shorten the operating life of the equipment, seriously damaging its components, especially the vaporizer, and invalidate the manufacturer's warranty.

11. STERILIZATION TABLE

Sharp surgical instruments must be wrapped to ensure they are sterile at time of use. The data shown in tables is only approximate: the type of sterilization cycle selected should be based on data supplied by the manufacturer of the items to be sterilized. This sterilizer is not suitable for the sterilization of liquids. The overall cycle time can change due to various factors (e.g. weight and type of load, etc.).

11.1 AUTOCLAVE B6 TYPE

Type and total duration of cycle	Cycle type	Steriliz. (Min.)	Drying (Min.)	Pressure	Max. load	Materials and instruments to be sterilised
Helix /B&D test 13 min	Test	3.30	2.30	2.06	Pack test only	Test cycle
Vacuum test 15 min	Test	/	/	-0.88	Nothing	Test cycle
121° Hollow wrapped 38 min	B	18.00	10.00	1.06	1 kg	Fragile hollow, stainless hollow and turbines (wrapped)
134° Hollow wrapped 24 min	B	5.00	10.00	2.06	1 kg	Stainless hollow and turbines (wrapped)
121° Hollow open 30 min	B	18.00	2.00	1.06	1 kg	Fragile hollow, stainless hollow and turbines (not wrapped)
134° Hollow open 21 min	B	5.00	2.00	2.06	1 kg	Stainless hollow and turbines (not wrapped)
PRION (134°) 44 min	B	20.00	10.00	2.06	1 kg	Stainless hollow and metal instruments (wrapped)
121° Porous 45 min	B	18.00	15.00	1.06	0.5 kg	Porous material (not wrapped)
134° Porous 31 min	B	5.00	15.00	2.06	0.5 kg	Porous loads (wrapped)
121° Solids open 31 min	S	18.00	2.00	1.06	1.8 kg	Rubber and fragile solids (not wrapped)
134° Solids open 21 min	S	5.00	2.00	2.06	1.8 kg	Rubber and metal solids (not wrapped)
134° Solids wrapped 26 min	S	5.00	10.00	2.06	1.8 kg	Rubber and fragile solids (wrapped)

THIS STERILIZER IS NOT DESIGNED FOR THE STERILISATION OF LIQUIDS

The times indicated in the timetable have been obtained trying to simulate very hard working conditions for the sterilizer. However, variations of the times indicated are possible: the supply voltage, the temperature at the beginning of the cycle, the weight of the load to sterilize and other conditions can have different effects on the total cycle time.

Note: cycle times indicated in the table have been obtained with warm sterilizer

11.2 AUTOCLAVE B18 TYPE

Type and total duration of cycle	Cycle type	Steriliz. (min.)	Drying (Min.)	Pressure	Max. load	Materials and instruments to be sterilized
Helix/B&D Test 22 min.	Test	3.30	4.00	2.06	Pack test only	Test cycle
Vacuum Test 15 min.	Test	/	/	-0.88	Empty	Test cycle
121°C Hollow Wrapped 45 min.	B	18.00	10.00	1.06	2.5 kg	Fragile hollow, stainless hollow and turbines (wrapped)
134°C Hollow Wrapped 34 min.	B	5.00	10.00	2.06	2.5 kg	Stainless hollow and turbines (wrapped)
121°C Wrapped Solids 42 min.	S	18.00	10.00	1.06	2.5 kg	Rubber and fragile solids (wrapped)
134°C Wrapped Solids 30 min.	S	5.00	10.00	2.06	2.5 kg	Rubber and metal solids (wrapped)
134°C Prion 49 min.	B	20.00	10.00	2.06	2.5 kg	Stainless hollow and metal instruments (wrapped)
121°C Porous 56 min.	B	18.00	15.00	1.06	1 kg	Porous material (not wrapped)
134°C Porous 46 min.	B	5.00	15.00	2.06	1 kg	Porous loads (wrapped)
121°C Fast 36 min.	S	18.00	4.00	1.06	4 kg	Rubber and fragile solids (not wrapped)
134°C Fast 25 min.	S	5.00	4.00	2.06	4 kg	Rubber and metal solids (not wrapped)
134°C Hollow Open 31 min	B	5.00	4.00	2.06	4 kg	Stainless hollow and turbines (not wrapped)

THIS STERILIZER IS NOT DESIGNED FOR THE STERILIZATION OF LIQUIDS

The times indicated in the timetable have been obtained trying to simulate very hard working conditions for the sterilizer. However, variations of the times indicated are possible: the supply voltage, the temperature at the beginning of the cycle, the weight of the load to sterilize and other conditions can have different effects on the total cycle time.

11.3 AUTOCLAVE B23 TYPE

Type and total duration of cycle	Cycle type	Steriliz. (min.)	Drying (Min.)	Pressure	Max. load	Materials and instruments to be sterilized
Helix/B&D Test 29 min.	Test	3.30	4.00	2.06	Pack test only	Test cycle
Vacuum Test 15 min.	Test	/	/	-0.88	Empty	Test cycle
121°C Hollow Wrapped 51 min.	B	18.00	15.00	1.06	3.5 kg	Fragile hollow, stainless hollow and turbines (wrapped)
134°C Hollow Wrapped 49 min.	B	5.00	15.00	2.06	3.5 kg	Stainless hollow and turbines (wrapped)
121°C Wrapped Solids 50 min.	S	18.00	15.00	1.06	3.5 kg	Rubber and fragile solids (wrapped)
134°C Wrapped Solids 40 min.	S	5.00	15.00	2.06	3.5 kg	Rubber and metal solids (wrapped)
134°C Prion 60 min.	B	20.00	15.00	2.06	3.5 kg	Stainless hollow and metal instruments (wrapped)
121°C Porous 65 min.	B	18.00	20.00	1.06	1.5 kg	Porous material (not wrapped)
134°C Porous 58 min.	B	5.00	20.00	2.06	1.5 kg	Porous loads (wrapped)
121°C Fast 44 min.	S	18.00	9.00	1.06	4.5 kg	Rubber and fragile solids (not wrapped)
134°C Fast 34 min.	S	5.00	9.00	2.06	4.5 kg	Rubber and metal solids (not wrapped)
134°C Hollow Open 41 min.	B	5.00	9.00	2.06	4.5 kg	Stainless hollow and turbines (not wrapped)

THIS STERILIZER IS NOT DESIGNED FOR THE STERILIZATION OF LIQUIDS

The times indicated in the timetable have been obtained trying to simulate very hard working conditions for the sterilizer. However, variations of the times indicated are possible: the supply voltage, the temperature at the beginning of the cycle, the weight of the load to sterilize and other conditions can have different effects on the total cycle time.

11.4 NIGHT CYCLE

If the sterilizer is not in use, it will switch to energy-saving mode with only the display background lit. Pressing any key (except the START/STOP button), the display will show the result of the last operation carried out (e.g. "END CYCLE" and "OPEN THE DOOR"). Any cycle can be made a "night cycle".

In the class B and class S autoclaves, during the night (and all the times a cycle ended), the vacuum pump remove the residual steam from the chamber every 30 minutes. This automatic operation serves to minimize the lasting condensate to better dry the instruments.

ATTENTION:

1. After completing the "night cycle" then opening of the door, it is normal to find condensed water on the gasket of the door and the front of the sterilization chamber.
2. In case of alarms (see paragraph 14) repeat the operation!

12. TEST CYCLES


To establish the proper check procedures and frequency, refer to the instructions issued by the relevant health authority.

12.1 "BOWIE & DICK" CYCLE (for B & S class autoclaves only)


The Bowie & Dick cycle verifies that the steam penetrates properly inside a porous load. Complete test cycle after removing all trays from the sterilizer chamber except the central one: place the test pack without any other instruments on this tray. Select the "Helix / B&D" cycle and start the cycle. The cycle result is confirmed by the test pack result.

12.2 "HELIX TEST" CYCLE (for B class autoclaves only)

This test should be performed to check steam penetration into a hollow load. Carry out the test cycle after removing all trays from the sterilizer chamber except the central one: place the "Helix-Indicator System" on this tray (without any other instruments). Select the "Helix / B&D" cycle and start the cycle. The cycle result is confirmed by the test paper included in the test material (see drawing below). **Important: the cycle must be made when the sterilizer is warm (immediately after having performed a work cycle).**

Bowie & Dick test

Helix test

12.3 VACUUM TEST CYCLE (for B & S class autoclaves only)

The Vacuum Test cycle is used to detect any pressure losses in the sterilization chamber. This test should be carried out on an empty machine, prior to starting other sterilization cycles.

The cycle cannot start if the temperature inside the sterilization chamber is $>40^{\circ}\text{C}$.

- Select "Vacuum Test" and start the cycle. The sterilizer reaches the selected vacuum degree and will hold it for 15 minutes.
- The test result is given by the "END CYCLE" and "OPEN THE DOOR" messages appearing on screen and by the relevant printout.
- If the machine failed the test (AL0600 or AL0601), the door seal should be checked, cleaned or, if necessary, replaced (paragraph 15.1).
- Similarly, the sterilization chamber rim should also be checked. Finally, the test should be repeated.

A failed Vacuum Test does not prevent the use of the sterilizer in the near future; however, it is recommended that you contact your technical support centre as sterilization cycles may be affected in the long term.

12.4 BIOLOGICAL INDICATOR TEST (for ALL autoclaves)

Together with the other chemical indicator tests, a biological indicator test may be required. This test involves sterilizing one or more vials containing biological spores, together with the normal items to be sterilized. On completion of the cycle, remove the vials and leave them to cool for a few minutes (following the manufacturer's instructions for the control procedure). Sterilized vials should normally be broken using the special tools supplied by the manufacturer and then inserted into a specific incubator: together with these, another vial that has not undergone sterilization (the Control) should be added for comparison. After the incubation period, the difference in colour in the sterilized vials will indicate whether the cycle was successful.

13. RECOMMENDATIONS FOR STERILIZATION

In order to prolong the life of both sterilizer components and instruments, recommended procedures should be followed as well as local health authority instructions. Below is a list of some precautions to be followed.

ATTENTION:

Do not use too much lubricating oil because it can cause damage to the valves of the vacuum pump and solenoid valves. This damage will not be considered under warranty.

1. Instruments should be cleaned with appropriate disinfectants immediately after use.
2. Brush the instruments to remove any residue.
3. Rinse the instruments in running water at room temperature.
4. Place the instruments in an ultrasonic bath treatment.
5. Rinse the instruments in demineralised water at room temperature. Thoroughly rinse and clean tools to eliminate possible residues of chemicals or lubrication oils that may damage components of the sterilizer.
6. Dry the instruments thoroughly.
7. Place the instruments on the sterilizer tray, so that the bags do not overlap. Use the tray holders and trays provided. Do not place tools in direct contact with the chamber.
8. If you need to sterilise instruments which are not wrapped in bags, the tray should be covered with suitable napkins, to ensure that each sterilised instrument is perfectly dried.
9. Instruments such as scissors or forceps should be opened slightly. Mirrors should be positioned face down.
10. Place bags with the paper side up.
11. If empty containers are sterilised, place them upside down to prevent water from accumulating.
12. At the end of the sterilization cycle, use the removal key provided as trays and tray holders are very hot.

The instructions listed above show how important the correct preparation of instruments is to successful sterilization. Even a single instrument carrying traces of disinfectant placed in the sterilizer may cause damage to the sterilization chamber and the instruments inside. Where, the sterilization process might be adversely affected even if no alarm is triggered.

Always use the appropriate indicator tests to verify successful sterilization of the instruments.


14. ALARMS AND ERRORS

Alarms visible on the display (paragraph 14.1) will prevent any subsequent operation: **in order to reset an alarm press SELECT and START buttons at the same time** until the display temporarily disappears. Alarms are also logged on the printout (see table below). Conversely, errors (paragraph 14.2) do not allow the cycle to restart, but indicate that a correct action is required before continuing with sterilization (e.g.: "DRAIN WATER").


ATTENTION:

Any non-completed cycle means that the instruments are not sterilized: on the display there is an alarm code (see the following pictures) and instruction for reset.

In the event of alarms, the cycle should be considered invalid and the load not sterilized

Interpreting ALARM codes:

14.1 ALARMS

Code and meaning	Problem	Solutions
		Reset alarm = SELECT+START at least 5 sec.
AL0001 Cycle stopped by user.	Occurs if the START/STOP button is pressed for one second or longer	Reset the alarm and repeat the cycle.
AL0002 No mains power supply	No mains supply: it is caused by the lack of main power supply.	
AL0003 Door open during cycle	Occurs if a control micro-switch detects an open door while the cycle is running.	Reset the alarm and repeat the cycle. If the problem persists, contact technical support.
AL0004 Timer stopped	The timer will stop when PCB battery power is insufficient.	Reset the time and date again following the instructions in paragraphs 8.2 and 8.3. Leave the sterilizer ON for at least one hour.
AL0005 High voltage	It is caused by an overload of the mains supply.	Reset the alarm and repeat the cycle.
AL0011 No 1st vacuum	The alarm is triggered only if the first vacuum is not reached.	Reset the alarm, clean the gasket then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact technical support.
AL0012 No 2nd vacuum	The alarm is triggered only if the second vacuum is not reached.	
AL0013 No 3rd vacuum	The alarm is triggered only if the third vacuum is not reached.	
AL0015 No vacuum when drying	Occurs if the pre-set vacuum degree is not reached during drying. <u>If this alarm is triggered, the machine has already completed the sterilization.</u>	Reset the alarm and repeat the cycle. Check if the problem persists in the following cycles: in this case contact technical support.
AL0021 No 1st pressure rise	The machine is unable to reach the 1 st pressure value set.	Reset the alarm, fill the tank until the maximum icon appears on screen and repeat the cycle: if the problem persists in the following cycles, contact technical support.
AL0022 No 2nd pressure rise	The machine is unable to reach the 2 nd pressure value set.	
AL0024 No final pressure rise	The machine is unable to reach the operating pressure.	
AL0031 No 1st pressure release	After reaching the first set pressure, the machine triggers an alarm.	Reset the alarm, remove the tray holder, clean the inside of the sterilization chamber and repeat the cycle.
AL0032 No 2nd pressure release	After reaching the second set pressure, the machine triggers an alarm.	
AL0034 No final pressure release	The machine fails to release pressure while drying.	
AL0100 T1 probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0101 OPEN T1	System detects that T1 probe is open.	
AL0102 T1 short circuit	System detects short-circuit in T1 probe.	
AL0110 High temp. in T1 probe	T1 probe has exceeded the pre-set cycle temperature.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact technical support.
AL0111 Low temp. in T1 probe during sterilization	Temperature detected by T1 probe during sterilization is below the minimum set value.	

Code and meaning	Problem	Solutions Reset alarm = SELECT+START at least 5 sec.
AL0200 T2 probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0201 OPEN T2	System detects that T2 probe is open.	
AL0202 T2 short circuit	System detects short-circuit in T2 probe.	
AL0210 High temp. in T2 probe	T2 probe has exceeded the pre-set cycle temperature.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact technical support.
AL0211 Low temp. in T2 probe during sterilization	T2 probe detects excessively low values during sterilization.	
AL0300 P probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0301 OPEN P	System detects that P probe is open.	
AL0302 P short circuit	System detects short-circuit in P probe.	
AL0310 High pressure during sterilization	P probe detects excessively high values during sterilization.	Reset the alarm, remove the tray holder, clean the sterilization chamber inside and repeat the cycle.
AL0311 Low pressure during sterilization	P probe detects excessively low values during sterilization.	Reset the alarm, clean the gasket and repeat the cycle: if the problem continues call technical support.
AL0400 TVP probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0401 OPEN TVP	System detects that TVP probe is open.	
AL0402 TVP short circuit	System detects short-circuit in TVP probe.	
AL0404 Low temp. in TVP probe	TVP probe does not reach operating temperature.	Reset the alarm and repeat the cycle. If the problem persists, contact technical support.
AL0405 High temp. in TVP probe	TVP probe has reached a temperature higher than operating threshold.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact technical support.
AL0500 TRS probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0501 OPEN TRS	System detects that TRS probe is open.	
AL0502 TRS short circuit	System detects short-circuit in TRS probe.	
AL0504 Low temp. in TRS probe	Temperature does not reach operating temperature.	
AL0505 High temp. in TRS probe	TRS probe has reached a temperature higher than operating threshold.	

Code and meaning	Problem	Solutions
		Reset alarm = SELECT+START at least 5 sec.
AL0600 Excessive drop during the VACUUM TEST stabilization phase	An excessive pressure drop occurred during the first 5 minutes of the stabilization phase, while performing the VACUUM TEST cycle.	Take the tray holder out: clean and <u>dry the sterilization chamber thoroughly</u> . Repeat the cycle: if the problem continues run a Vacuum Test cycle when the sterilizer is cold (see paragraph 12).
AL0601 Excessive drop during the VACUUM TEST maintenance phase	An excessive pressure drop occurred during the first 10 minutes of the maintenance phase, while performing the VACUUM TEST cycle.	
AL0700 T1/T2 comparison	Conflicting temperature values detected by the two interior sensors during sterilization.	Reset the alarm, fill the tank until the maximum icon appears on the screen and repeat the cycle: if the problem persists in the following cycles, contact technical support.

14.2 ERRORS

The following table contains a list of messages which may appear on display:

MESSAGE	CAUSE	SOLUTION
OPEN THE DOOR	At the time of power-up you are prompted to open the door so that the machine can perform a pressure check.	Open the door to allow automatic setting of atmospheric pressure (two beeps).
DOOR OPEN	A cycle has been started with the door open.	Close the door and start the cycle.
DOOR NOT LOCKED	Even though the door is perfectly closed, it is read as "not locked" due to mechanical reasons.	Press the DOOR button, open and close the door, then press the START button to start the cycle. If this error persists, contact technical support
TO LOAD WATER	A cycle has been started with the clean water minimum level message flashing on the display.	Fill the tank with clean water (paragraph 10.1).
TO DRAIN WATER	A cycle has been started with the dirty water maximum level message flashing on the display.	Drain dirty water from the tank (paragraph 10.4).
FULL WATER RESERVOIR	Attempt to fill the tank with clean water with the full tank icon already present on the display.	Disconnect the water hose, then start a new cycle.
DOOR STILL LOCKED	Door lock system fails to open completely when cycle terminates.	Close the door again and start a new cycle: stop the cycle after a few seconds with the START/STOP button (paragraph 6.2), reset the alarm (paragraph 14), then try opening the door with the DOOR button. If the error persists, contact technical support.
CHAMBER TEMP.>40°	Attempt to start a VACUUM TEST cycle with a chamber temperature higher than 40°C; wait for temperatures to drop before starting the test cycle.	Open the door and wait for the temperature reading on the display to fall below 40°C.
"POOR WATER" or "BAD WATER"	OPTIONAL: the autoclave detects that water quality is just sufficient (POOR) or scarce (BAD).	See chaps. 10.1, 10.2 and 10.6.
BACT. FILTER	The sterilizer cycle counter prompts the operator to <u>change</u> the bacterial filter.	Replace the used bacterial filter with a new one (see paragraph 15.3). Once replaced, with the door open, reset the message pressing START for 10 seconds, at least, until the long beep.
SERVICE 1xx (& 2xx & 3xx)	The machine has performed a high number of cycles and requires the intervention of a technician for a periodic inspection.	Contact your dealer as soon as possible, indicating the exact code provided by on the machine and requesting technical assistance.

15. MAINTENANCE

Disconnect the equipment from the mains before carrying out any maintenance.
Make sure that the sterilizer has cooled down before accessing the chamber or parts inside.

15.1 DAILY MAINTENANCE

Daily maintenance includes keeping the door seals in good working order, cleaning the chamber rim area (very important for a successful test cycle) and checking the water levels in the tanks.

- **DOOR SEAL:** Clean the door seal using the soft part of the sponge supplied with the equipment. Cleaning should be carried out to remove any impurities which might affect the test cycles.
- **CHAMBER RIM:** It is the outer rim of the sterilization chamber with which the seal makes contact. Use the coarse part of the sponge supplied.
- **WATER LEVELS** (paragraph 6.3): Before starting another sterilization cycle, check the water level in the tanks.
- **GENERAL SURFACE CLEANING:** Use a cloth to remove dust and other deposits from the top of the machine.

15.2 WEEKLY MAINTENANCE


Weekly maintenance requires a visual check and cleaning of the inside of the sterilization chamber. Remove the trays and tray holder from the chamber before cleaning.

- **CHAMBER INTERIOR:** Use the coarse part of the sponge supplied to remove small impurities from the bottom of the chamber. If lime deposits are found, it is advisable to check the quality of water being used. Do not use disinfectants or sharp objects.

15.3 QUARTERLY MAINTENANCE

Quarterly maintenance requires lubrication of the door hinges and replacement of the bacterial filter.

- **HINGE LUBRICATION:** Spray the two door hinges with small quantities of silicone oil.
- **REPLACE THE BACTERIAL FILTER:** Replace the filter when required by the autoclave and reset the counter (see paragraph 14.2). A clogged bacterial filter may lead to difficulties opening the door, on completion of the sterilization cycle and reduces the sterilizer drying performance.

15.4. ANNUAL MAINTENANCE/VALIDATION

The sterilizer is fundamental to the protection of the patient and operator alike. Even though the electronic controls of these machines are increasingly reliable, it is good practice to carry out a functional check of the equipment at least once a year. This check should be carried out by authorised service centres only, with calibrated and certified instruments, in order to guarantee the long life and reliable operation (validation) of the equipment. To establish the proper check procedures, refer to the instructions issued by the relevant health authority.

- **ANNUAL VALIDATION:** Requires the use of instruments calibrated by specialist centres to check the sterilization cycle parameters. This check includes inspections of pressure and temperature probes and the timer. On request, the manufacturer issues a yearly test certificate for machines returned to its premises for maintenance and checks.
- **PRESSURE VESSEL MAINTENANCE:** in addition to the annual inspection indicated by the sterilizer via its "service light," checks on the following must also be carried out: the closure system, safety valve, and seals. As this is a pressure vessel, these checks must be carried out by specialized technicians using original spare parts. Furthermore, perform a visual inspection on the pressure vessel at least every 2000 cycles (corresponding to approximately 20% of the useful life of the device, evaluated as 10,000 cycles) to verify the absence of cracks or other critical issues that compromise safety.

16. TECHNICAL FEATURES

MECHANICAL SPECIFICATIONS

Operating temperature	+5° +30°C		
Max. operating altitude	2.000 m		
Max. relative moisture at 30°C	80%		
Max. relative moisture at 40°C	50%		
Overall dimension:	6 Lts.	18 Lts.	23 Lts.
B type	450 x 345 x 610 mm	505 x 400 x 615 mm	505 x 400 x 690 mm
S & N type	450 x 345 x 610 mm	445 x 400 x 615 mm	445 x 400 x 690 mm
Weight (empty)	6 Lts.	18 Lts.	23 Lts.
Type/kg	B & S=35	B=45 / S=40 / N=35	B=50 / S=45 / N=40
Weight (full tanks)	6 Lts.	15 & 18 Lts.	23 Lts.
Type/kg	B & S=40	B=50 / S=45 / N=40	B=55 / S=50 / N=45
Overall dimensions (door open)	6 Lts.= 310 mm 18 & 23 Lts.= 350 mm		
Weight per unit area of support	Max 32,46 N/cm ²		
Volume	Max 0.14 m ²		
Colour	RAL9002		
Material	AISI 304 / FeP01		
Noise level at 1 mt distance	53,6 dbA		
Noise level in front of the display	62,2 dbA		

ELECTRICAL SPECIFICATIONS

Supply voltage	230 Vac (+/-10%)		
Power (Watts)	B6 Lts. = 2500	S6 Lts. = 1500	B/S/N 18/23 Lts. = 2400
Frequency	50-60 Hz		
Supply cable (L 1.5m)	2+1 x 1,5 mm ²		
Fuses	6,3x32mm - T12A		
Max. transmitted heat	5,76 MJ/h (1370 Kcal/h)		
Insulation class	1		

CHAMBER SPECIFICATIONS

Max. operating pressure	2.5 bars (relative)		
Max. operating vacuum	-0.90 bars (relative)		
Max. operating temperature	140°C		
Material	Stainless steel AISI 304		
Dimensions: ØxP	6 Lts 170 x 265 mm	18 Lts. 236 x 381,5 mm	23 Lts. 236 x 530 mm

CLEAN WATER TANK SPECIFICATIONS

Max. capacity	4,5 Lts.		
Cycles per tankful	6 Lts.= ~12	18 Lts.= ~6	23 Lts.= ~5
Material	Polyethylene		

DIRTY WATER TANK SPECIFICATIONS

Max. capacity	4,5 Lts.		
Cycles per tankful	6 Lts.= ~12	18 Lts.= ~6	23 Lts.= ~5
Material	Polyethylene		
Maximum temperature of drained water	50°C		

BACTERIAL FILTER DATA

Maximum diameter	56 mm		
Max. allowed filtering capacity	0.3 micron		
No. of cycles before replacement	~300		

TRAY HOLDER SPECIFICATIONS

Material	Stainless steel		
----------	-----------------	--	--

TRAY SPECIFICATIONS

Material	Anodised aluminium		
----------	--------------------	--	--

CLASS OF POLLUTION	2		
--------------------	---	--	--

THE MANUFACTURER RESERVES THE RIGHT TO IMPLEMENT TECHNICAL MODIFICATIONS WITHOUT PRIOR NOTICE. THIS MANUAL IS THE SOLE PROPERTY OF MANUFACTURER: BY LAW, IT CANNOT BE REPRODUCED OR TRANSFERRED TO THIRD PARTIES WITHOUT ANY WRITTEN AUTHORISATION.

17. WARRANTY

For all defects of conformity existing at the time of delivery of the unit and attributable to actions or omissions by the manufacturer

A. The manufacturer guarantees this product for a period of:

1. 24 months or 1500 cycles on the whole product and
2. 5 years or 5000 cycles on the chamber only.

In both cases (1 and 2) the warranty expires when one of the two aforementioned conditions is met. Labour is excluded from the warranty, if it does not take place at the manufacturer's site.

The warranty period begins on the date the machine is delivered to the customer. In the event of any contestation, the date of delivery is considered valid as evidenced by a document valid for tax purposes (packing slip, invoice, tax receipt or similar) which must contain the name of the seller, the date of delivery, the identification details of the product (serial number and model) and the sales price. The installation documents present inside the product must be completed in full, stamped and signed by both the dealer and the customer and sent to the manufacturer, otherwise the warranty will not be effective.

B. In order for this WARRANTY to be fully valid, the following is necessary:

1. All installation operations must have been carried out strictly following the instructions provided in this USER INSTRUCTION MANUAL.
2. All use and routine maintenance must be carried out according to the USER INSTRUCTIONS.
3. All planned "SERVICE" interventions that the product needs and indicates must be carried out. The interventions needed to solve the "SERVICE" procedures are always excluded from the warranty.
4. Any repair carried out under the warranty must be performed by authorised personnel and only original spare parts used. Parts replaced under the warranty must be returned to the supplier (or they will be charged for), with the exception of agreements stipulated in advance between the parties.

C. THE FOLLOWING ARE NOT COVERED BY THE WARRANTY:

1. Labour, when this does not take place at the manufacturer's site.
2. Damage resulting from transportation, with the exception of agreements stipulated in advance between the parties.
3. All components that manifest a lack of conformity resulting from incorrect installation of machinery.
4. Damage caused by poor maintenance, neglect or carelessness of usage by the user and failure to comply with what is set out and recommended in the USER INSTRUCTIONS booklet.
5. Damage caused by tampering with the product or product parts.
6. Damage resulting from all other causes not attributable to the manufacturer.
7. All components subject to normal wear and tear (e.g. polycarbonate keypad, tubes supplied, seals, trays, filters etc.) and other accessories if it cannot be proven that it is the result of a manufacturing defect.
8. Costs of delivering spare parts and/or finished products.

D. LIMITATIONS OF THE WARRANTY:

1. The purchaser is not entitled to the replacement of the complete machine if the defect is not reported within two months of the date of purchase.
2. It is at the discretion of the manufacturer whether to repair or replace a component under the warranty. This does not include, however, the cost of labour and travelling expenses for personnel.
3. No compensation is awarded for machine downtime.
4. The warranty is automatically considered null and void if the machine is tampered with, repaired or modified by purchaser or third parties not authorised by the manufacturer. For interventions, the purchaser must contact the dealer or service personnel from the manufacturer only.

The manufacturer disclaims all liability for any damage or injury caused, directly and/or indirectly, to persons, property and animals due to failure to observe the general safety conditions and requirements set out in the USER INSTRUCTIONS MANUAL, in particular concerning the instructions for product installation, use and maintenance.

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BMS DENTAL Srl